

	Salt Oven Usage Procedures	
Document Number:		Revision Number: 01
Document Owner: SNO+ Calibration Post-Doc		
Reviewer:		
Name:	Signature	Date
Authorizer:		
Name:	Signature	Date

1 Purpose

This document describes the procedure for preparing and operating the Barnstead Thermolyne F21100 for use with a Sodium Chloride (NaCl) loaded quartz tube. The tube furnace heats the salt to a desired temperature of 610C, just below its melting point. A stream of Helium flows through the heated tube, along with the heated NaCl aerosol. The tube has a water cooling section before the gas outlet, preventing damage to the following gas lines and the Umbilical.

Warnings Refer to THA for required precautions needed for tube furnace operations. Always use heat resistance gloves when handling the salt tube. Ensure a proper barricade is in place around the oven work space. Ensure proper signage is in place to warn others of the hazards present beyond the barricade.

2 Preparing Oven

- 1: Only use Tube Furnace in designated area
- 2: Make sure the furnace is empty before loading the salt tube
- 3: Load the salt tube in the oven
- 4: Connect Helium gas inlet and outlet to the tube
- 5: Connect the output end to the water cooling loop using the attached flex hose

3 Starting Water Cooling Loop

- 6: Ensure the needle valve behind the gas board is closed.
- 7: Check that the inlet line is attached to the water line after the valve and that the outlet line is running into the drain channel.
- 8: Insert the input line into the sleeve on the salt tube closer to the oven. Clamp the line using the hose clamp provided.
- 9: Insert the output line into the sleeve on the salt oven nearer to the He gas outlet. Clamp the line using the hose clamp provided.

Salt Oven Usage Procedures	Document No:
	Revision No: 01
	Effective Date: 2018-04-30
	Page 2

10: Slowly open the valve. Ensure that all of the bubbles are out of the line before starting the Helium flow.

11: Check the seal on the input and output lines for water leaks. If the water pressure is too high water will leak from the cooling loop connections. If there is a leak close the water valve until the leak is no longer present.

4 Turn on Oven

12: Turn on the furnace, make sure the setpoint is 615C

13: Use the “UP” or “DOWN” button while in the home display to change the setpoint

4.1 Tuning Cycle

14: Press the “PAGE” button until “Atun” is displayed

15: Press the “SCROLL” button, the display will read **tunE**

16: Press either “UP” or “DOWN” to select **ON**

17: Press both “SCROLL” and “PAGE” to return to the main display

18: Wait for the tuning to complete

4.2 Continuous Heating

19: Press the “PAGE” button until **SP** is displayed

20: Press the “SCROLL” button until **tm.OP** is displayed

21: Press either “UP” or “DOWN” to select **Opt.2**. This setting will ramp and maintain setpoint temperature

22: Press the “SCROLL” button until **StAt** is displayed

23: Press either “UP” or “DOWN” to select **run**

24: Press both “SCROLL” and “PAGE” to return to the main display

Salt Oven Usage Procedures	Document No:
	Revision No: 01
	Effective Date: 2018-04-30
	Page 3

5 Flow Gas

- 25: Ensure that all of the valves on the gas board are closed.
- 26: Check that the valves on the He, N₂, and CO₂ bottles are closed.
- 27: Verify that the "Solenoid Valve II" switch is in the OFF position.
- 28: Enter the DCR using standard entry procedure and verify that
 - the gas input line is connected to the "dry" end of the umbilical
 - the gas return line is connected to the "dry" end of the umbilical.
- 29: Turn on the power onf the main instrument panel and the dual flow controller box.
- 30: Flip the switch on the flow controller box to "B" (for N₂ and/or CO₂). The reading takes a few minutes to equilibrate to zero or near zero.
- 31: Ensure that the CO₂ and He bottles are hooked up to the input line at VA2. Ensure that the N₂ bottle is connected to VA1.
- 32: Check that VA2 and VA2 are closed.
- 33: Close the needle valve on the CO₂ bottle.
- 34: Open the main valves on the CO₂, He and N₂ bottles.
- 35: Set the regulator on the CO₂ and He bottle to 33 psi. Set the N₂ regulator to 20 psi.
- 36: Record the bottle pressure:

Transducer	Reading
CO ₂ Bottle Pressure	
He Bottle Pressure	
N ₂ Bottle Pressure	

The pressure of a full CO₂ bottle is 850 to 900 psi. The CO₂ in the bottle is in liquid or solid form and the pressure will stay relatively high until gas only and will drop rapidly thereafter. The helium bottle is simply gas under pressure and the pressure will drain steadily throughout its use. ...

6 Turn off Oven

- 37: Press the "PAGE" button until **SP** is displayed
- 38: Press the "SCROLL" button until **StAt** is displayed

Salt Oven Usage Procedures	Document No:
	Revision No: 01
	Effective Date: 2018-04-30
	Page 4

39: Press either “UP” or “DOWN” to select **off**

40: Press both “SCROLL” and “PAGE” to return to the main display

41: The oven should start to lose heat, let the salt tube cool to ambient temperature

42: Turn off tube furnace